

# WILD M3B · M3C · M3Z · M3Z PLAN

Instructions for use  
Mode d'emploi  
Bedienungsanleitung  
Instrucciones para el uso

 **WILD LEITZ**

## Instructions for use

### Instruction checklist

#### I. Operation

1. Interpupillary distance
2. Exit pupil position
3. Magnification changer
4. Focusing
5. Dioptic correction
6. Drive housing type "S"

#### II. Assembly

7. Stands
  - 7.1 Swinging-arm and table-clamp stands
  - 7.2 Swinging-arm stand, light
  - 7.3 Universal stand
8. Binocular tubes
9. Objectives
10. Eyepieces

#### III. Illuminators

11. Incident light
  - 11.1 Mains lamp 25 W
  - 11.2 Low-voltage lamp 6 V/10 W
  - 11.3 Low-voltage lamp 6 V/20 W
  - 11.4 Ring illuminator
  - 11.5 Near-vertical illuminator
  - 11.6 Coaxial incident illuminator
12. Transmitted light
  - 12.1 Transmitted-light stand, bright field, circular
  - 12.2 Transmitted-light stand EB
  - 12.3 Transmitted-light stand, bright/dark field
13. Combined illuminators
14. Regulating transformer

The stereomicroscope produces a spatial image of three-dimensional objects. Two separate beam paths which traverse a common main objective cause each eye to observe the object from a different direction. The two images are fused into a single plastic image.

Wild stereomicroscopes are designed with two separate optical trains, the axes of which run parallel to one another after passing through a common main objective. Because the intermediate image planes are parallel to the object plane, the whole field is uniformly sharp and the natural link between accommodation and convergence is retained.

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#### **IV. Accessories**

- 15. Filterholder
- 16. Polarisation
  - 16.1 Polariser on glass stage plate
  - 16.2 Rotatable Pol. stage
- 17. Stages
  - 17.1 Gliding stage
  - 17.2 Cup stage
  - 17.3 Stage carrier with mechanical stage
  - 17.4 Transmitted-light stage
- 18. Handrests
- 19. Attachment for vertical and oblique observation
- 20. Measuring
- 21. Tubes
  - 21.1 Double-iris diaphragm
  - 21.2 Discussion tube
  - 21.3 Drawing tube
  - 21.4 Phototubes
  - 21.5 Beam splitter

#### **V. Appendix**

- 22. Choosing the optics combination
  - 22.1 Working distance
  - 22.2 Total magnification
  - 22.3 Field diameter
  - 22.4 Resolution
  - 22.5 Depth of field
- 23. Care of the instrument
- 24. Dimensions
- 25. Optical data

**The multiple-coated optics of a WILD Stereomicroscope result in first-class image transfer characteristics and consequently in a crisp image.**

**The human-engineered design ensures a comfortable and fatigue-free working position with conveniently-situated controls. A rich choice of modular accessories solves many problems.**

**The world-renowned Swiss quality of the Wild instruments guarantees the long-term untroubled functioning of your stereomicroscope and is supported by the technical advice and the servicing facilities available in the Wild agencies which are established in over one hundred countries.**

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
### Here's how to use your Wild Stereomicroscope properly:

1. Set up the stand.  
Slide your stereomicroscope down the column. Sect. 7
2. Fit the binocular tube and additional tubes if required. Sect. 8  
Sect. 21
3. Insert the eyepieces.  
Secure them with the clamping screws. Sect. 10
4. Choose an objective.  
Set the working distance. Sect. 9
5. Mount the illuminator required and switch it on. Sects. 11-14
6. Look down the binocular tube. Displace the eyepiece tubes to give the correct interpupillary distance. Sect. 1
7. Set the dioptic correction for both eyes by turning the diopter rings of the eyepieces. Sect. 5
8. Position the object to be investigated under the objective.  
Use the focusing drive to bring the image into focus. Sect. 4
9. Select a low magnification in order to obtain an overall view of the object.  
Rotate the magnification changer to the desired magnification. Sect. 3

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




# I. Operation

## 1. Interpupillary distance

The interpupillary distance can be individually set on all binocular tubes by displacing the two tubes within the range from 52 mm to 76 mm (fig. 1).

If the interpupillary distance has been correctly set, a single circular field is seen, which results from the fusion of the two separate fields in the stereomicroscope.

## 2. Exit pupil position

The correct distance of the eye from the eyepiece is important in the quality of the image. The further the eye is from the correct position (the eyepoint), the greater the loss of the peripheral part of the field, i.e. the field diameter is reduced and the marginal areas become dark.

At the eyepoint, all of the ray bundles emerging from the eyepiece are combined (fig. 2).

This position is located by moving the eye slowly towards the eyepiece until the whole field of view is seen sharply defined.

The stereomicroscopes of the M3 series are routinely supplied with two 10×/21B distortion-free wide-field eyepieces for spectacle wearers (fig. 3). Their exit pupil distance of 22 mm permits work either with or without spectacles.

The degree of eye contact can be adjusted by displacing the eyecups. Keep hold of the diopter ring when pulling out the eyecups.

To extend the magnification range, adjustable wide-field eyepieces with intrinsic magnifications 15×, 20× and 32× are available; they have a 12 mm exit pupil distance and are provided with attachable eyecups (fig. 4), but are not suitable for spectacle wearers, who should instead use the 15× and 20× eyepieces for spectacle wearers (fig. 5). These latter eyepieces have a 22 mm exit pupil distance and screw-on eyecups. Users who prefer closer eye contact can fit 6 mm or 9 mm spacing rings between eyepiece and eyecup.

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
### 3. Magnification changer

When the magnification changer (2) is turned, the magnification is altered by the following factors:

- M3B by 0.64, 1.6 and 4.0
- M3C by 0.64, 1.0, 1.6, 2.5 and 4.0
- M3Z steplessly within the factor range 0.65 to 4.0.

The focusing knobs are engraved with the total magnifications obtained, assuming that a 1.0× objective and 10× eyepieces are used. Additional optical data is given in the table (page 78).

Click-stops for the total magnifications 10×, 16× and 25× are engaged when the lever (6) on the model M3Z is pointing towards the observer.

### 4. Focusing


The bilateral focusing knobs (1) ensure comfortable operation. The whole length of the side-faced column is available for the focusing procedure.

#### Adjust ease of movement of focusing drive:

- Rotate the drive knobs until the concealed screws at the rear are visible through the cutouts.
- Use an Allen key to change the positions of the hollow screws. They should both be turned by about the same angle.

The **inclinable drive housing** has a focusing range of 50 mm (fig. 3).

The ease of movement is individually adjustable by turning the two drive knobs in opposite directions.

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
## 5. Dioptric correction

With the diopter rings of the eyepieces, the observer can correct long or short sightedness over the range  $\pm 5$  diopters.

If the correction has been performed properly, the image will stay sharp for both eyes when the magnification is changed.

### Preparations

The drive housing "S" is to be positioned for stereoscopic viewing (see section 6).

- a Set the working distance (see section 9).
- b Set the interpupillary distance (see section 1).
- c Set the diopter rings of the eyepieces to "0" (c).
- d Open the double-iris diaphragm (if present) (see section 21.1).

- Set the magnification changer to the highest position.

**Focus on a flat test specimen** (e.g. a cross on paper) **with the focusing drive (1).**

- Set the magnification changer to the lowest position.

Turn the diopter rings of the eyepieces fully out in the '+' direction without looking into the eyepieces.

**Turn the diopter rings of the eyepieces individually in the '-' direction** until the specimen appears sharp. While doing this, close the other eye.

- Set the magnification changer to the highest position.

**Adjust the sharpness** of focus and correct it with the focusing drive if necessary.

The specimen should now remain in focus at all zoom positions without any further refocusing. Otherwise the procedure must be repeated.

The diopter values set remain the same irrespective of the objective or eyepiece used.

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
## 6. Drive housing type "S"



Any optics carrier from the WILD M3 Series can be equipped either with the standard drive housing or with the version "S".

The Drive Housing Type "S" offers the choice between vertical viewing down either the left or the right beam path, or stereoscopic viewing using both beam paths simultaneously.

The stereo position is necessary for **observing three-dimensional objects** and for work on the object.

The vertical-beam position is to be preferred for **photography and for measuring**, and this configuration is also better for the study of polarisation colours.

- By pressing knob (7), the optics carrier can be moved into the required position.
- For photography and for measurement, the beam path to be chosen is the one above which the camera or the measuring eyepiece is positioned.
- The photo outlet of the inclined trinocular tube is over the right beam path. Move the optics carrier to this position (fig. 3).

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

## II. Assembly

The modular design of the M3 Series enables a stereomicroscope outfit to be assembled in accordance with customer requirements. The optics carrier and the drive housing are available in various versions and can be combined as desired.

### Optics carriers:

- a) M3B with three-step magnification changer
- b) M3C with five-step magnification changer
- c) M3Z with 1:6 zoom magnification changer

### Drive housings:

- d) Standard drive housing for stereoscopic observation
- e) Drive housing type "S" for stereoscopic and vertical observation
- f) Drive housing, inclinable, for swinging-arm and table-clamp stands and for OEM adaptations.

The complete outfit comprises the optics carrier, drive housing, binocular tube, objectives and eyepieces. These components are supplemented by the particular stand and illuminator chosen, and by any special modular accessories which may have been selected.

The M3B/C/Z Stereomicroscopes are delivered with achromatic objectives; the M3Z PLAN with a planachromatic objective.

The M3Z Stereomicroscope is delivered with the magnification changer engaged at the 16 $\times$  position. This position is to be engaged if the instrument is later transported (see section 3).

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


## 7. Stands

**Incident-light stands are used for examining and manipulating three-dimensional objects in conjunction with lighting from above.** (Appropriate illuminators are indicated in section 11.)

The standard outfit includes the trapezoidal incident-light base with 320 mm side-faced column (fig. 1).

The side-faced column is secured to the base by three hollow screws.

- Place the 100 mm diam. stage plate in the port provided in the baseplate. Depending on the characteristics of the specimen, use either the black or the white side.

If an 80 mm diam. stage plate is used, combine it with a spacing ring.

- Fit the stage clips.

The adapter (36) serves to attach the side-faced column to the rectangular incident-light base (fig. 2). This base has an 80 mm diameter stage plate.

**Transparent or translucent objects such as insects, cultures, textiles or prepared specimens are investigated in combination with one of the transmitted-light stands (fig. 3).**

The transmitted-light stands are described in section 12.

The side-faced column is fitted to the transmitted-light stands by means of an adapter:

- Secure the adapter (36) to the side-faced column by means of the three hollow screws.
- Connect the side-faced column and the adapter to the transmitted-light base by means of the two additional hollow screws.
- Fit the 80 mm diam. stage plate and the stage clips.

Slide the optics carrier/drive housing combination on to the side-faced column so that the rack and pinion engage and the focusing drive (1) can be operated. The side-faced column ensures that the instrument is always precisely aligned to the base.

A 445 mm long side-faced column is available for longer working distances. It is usable on incident- and transmitted-light stages (with the exception of the transmitted-light stand with circular base).

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

### 7.1 Swinging-arm and table-clamp stands

Swinging-arm stands and table-clamp stands are available for WILD stereomicroscopes of the M3 Series (figs. 1, 2). Each consists of a stable base or a clamp, a post (20) and a horizontal arm (21).

In the standard outfits, the inclinable drive housing forms part of the swinging-arm stand; the side-faced carrier rod forms part of the table-clamp stand.

The post and the horizontal arm of the large swinging-arm stand and of the large table-clamp stand are longer.

#### Fitting

The post of the **swinging-arm stand** is secured to the base with four hollow screws.

The clamp of the **table-clamp stand** can be secured to table tops 20 mm–50 mm thick.

- Secure the cross-member (29) to the post with the clamping lever (24) and tighten the safety ring (23).
- Slide the horizontal arm into the cross-member and secure it.

#### Swinging-arm stand

- Secure the optics carrier to the inclinable drive housing, using the two hollow screws.
- Screw the objective on to the drive housing.
- Peg the drive housing into the coupling of the horizontal arm.
- Release the clamping lever (18) and adjust the tilt as required.
- Attach the lamp to the drive housing.

#### Table-clamp stand

- Slide the side-faced carrier rod over the horizontal arm and secure it with the clamping screw (26).
- The clamping lever (27) enables the side-faced carrier rod to be fixed at any angle. The two clamping levers can be pulled out along their axes and reset at another position.
- Remove the end-stop (28).
- Using two hollow screws, secure the optics carrier to the drive housing.
- Slide the stereomicroscope with drive housing on to the side-faced column so that the focusing drive (1) can be used. Fit the lampholder if required. Replace the end-stop.
- If the two hollow screws are loosened, the side-faced carrier rod can be turned by 90° on the tiltable joint.

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
### 7.2 Swinging-arm stand, light

The socket (37) is secured to benches by four hollow screws.

Place the column in the socket. The complete stand, with optics carrier attached, can easily be lifted out again later and remounted at another working place. Fit the clamping block (29) to the column.

Fit the cross-member (38) to the clamping block and secure it with the lever (26).


- Secure the optics carrier to the objective adapter with swinging arm (39) by means of two hollow screws.
- Screw on an objective of long focal length (e.g. the 0.63 $\times$ ).
- Fit the mounted optics carrier with swinging arm into the cross-member and tighten the lever (27).

### 7.3 Universal stand

The horizontal arm (21) of the large swinging-arm stand and the inclinable side-faced carrier rod (22) can be attached to the universal stand (fig. 2). The 50 mm diameter column (20) is available with lengths 450 mm or 800 mm and is secured with three hollow screws to the 52 cm $\times$ 34 cm baseplate.

Fit the safety ring (23), the horizontal arm (21), the inclinable side-faced carrier rod (22) and the stereomicroscope (see also section 7.1, table-clamp stand).

Stage carriers with magnetic linkages (90) for mechanical stages (fig. 2) and for special stages with 80 mm diameter mounts (see section 17) can be attached to the baseplate.

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
## 8. Binocular tubes

The standard outfit includes the 45° **inclined binocular tube** (fig. 1, right).

For outfits with transmitted-light stands or additional tubes, the low **inclined binocular tube** is recommended (fig. 1, left).

The **straight binocular tube** serves to ensure a comfortable working position when the instrument is inclined (fig. 2).

The low inclined binocular tube and the straight tube are also available as **version T** with built-in double-iris diaphragm and with drive for adjusting interpupillary distance.

A **tube wedge** (fig. 1, right) enables the inclination of the binocular tubes to be increased or decreased by 15°.

The **binocular tube with variable angle of observation** provides for fatigue-free body positions with all types of instrument assembly and with any combination of accessories (fig. 3).

The angle of observation can be steplessly varied by up to 90° in either direction from the vertical and can also be set for either of two viewing heights. The image automatically remains erect and laterally correct. This tube is not usable on phototubes.

### Fitting the binocular tube:

- Ease it into position so that the locating screw engages the slot.
- Tighten the clamping screw (4).

**Accessory tubes**, such as the beam splitter or phototube, are fitted between the optics carrier and the binocular tube in the same way (see section 21).

## 9. Objectives



**Interchangeable achromatic objectives** are available for the M3B/C/Z Stereomicroscopes. They enable the working distance, total magnification and field diameter to be changed (see table on p. 78).

The objectives are colour-coded.

The WILD M3Z PLAN Stereomicroscope has an 1.0× or 1.6× **planachromatic objective**. The adapter ring can be unscrewed and removed to increase the working distance (fig. 4).

### Fitting:

- Screw the required objective on to the drive housing in an anticlockwise direction.
- Using the focusing drive, set the correct working distance.

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
## 10. Eyepieces

Stereomicroscopes of the WILD M3 series are routinely supplied with two distortion-free 10×/21B eyepieces for spectacle wearers (fig. 1). The adjustment of the eyecups is described in section 2.

The magnification range can be extended by using 15×, 20× and 32× wide-field eyepieces, with plug-on eyecups (10, fig. 2), and 15× and 20× eyepieces for spectacle wearers, with screw-on eyecups (11).


All types of eyepiece have clamping screws (14) for securing them in the tubes.

**Graticules in mounts** fit into the eyepieces; they are used for photomicrography and for measuring (fig. 3):

- Unscrew the sleeve (16) from beneath the eyepiece.
- Fit the graticule and replace the sleeve.

The use of the graticules is described in the appropriate instructions (for MPS systems or for measuring).

Please also observe the notes in section 22 about optics combinations.

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


#### Subject

#### Incident light

- ① **Inclined incident light**, dark field
- ② **Incident ring light**, daylight quality, free of shadows and reflections
- ③ **Coaxial incident light** (following the observation beam path)

Opaque, irregular, rough, highly diffusing, mostly unprepared, e. g. plants and insects  
 Contrast-rich, deeply-fissured, unprepared, e. g. rocks, metal objects

Flat, polished and highly-reflecting surfaces, e. g. integrated circuits, metal sections

#### Near-vertical light

(impinging steeply from above)

Depressions, holes, bores, cavities, e. g. printed-circuit boards, hybrid circuits, jets, filters, cracks

#### Transmitted light

- ④ **Transmitted light, bright field**, directly through the specimen
- ⑤ **Transmitted light, dark field**, obliquely through the specimen
- ⑥ **Combined incident and transmitted light**
- ⑦ **Polarised transmitted light**

Transparent, translucent, mostly prepared, e. g. thin sections, polished sections, smears  
 Low-contrast, e. g. anatomical and histological sections  
 Partly-transparent material, e. g. textiles, grids, insects  
 Double-refracting material, e. g. thin sections of rocks and bones, plastics, liquid crystals

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
## 11. Incident light

The mains (line) lamp and the low-voltage lamps are used for the **inclined incident illumination** (dark field) of spatial subjects (fig. 1). To avoid too-strong shadows, the subject is best illuminated from both sides.

### 11.1 Mains lamp 25W

The 25 W mains (line) lamp is suitable for observations with Wild stereomicroscopes; for more demanding requirements and for photomicrography the 6 V/10 W and 6 V/20 W low-voltage lamps are recommended. The 25 W mains lamp is supplied for voltages of 110–130 V (115 V reflector bulb) and of 220–235 V (220 V reflector bulb). Make sure that you have the correct combination. The lamp housing (50) bears a rotatable coupling to attach the lampholder.


The mains (line) lamp is generally fitted to the objective mount by means of the articulated arm (45) and the adapter (48), or to the cast foot by means of the adapter (49).

The clampable lampholder (46) and the adapter (53) together enable the lamp to be secured to the side-faced column (see section 11.2).

The coupling system between lamp and lampholder is the same as between lampholder and adapter.

### Fitting to the objective mount (fig. 2)

- Unscrew and remove objective.
- Slide adapter (48) and spacing ring over objective.
- Replace objective.
- Secure lampholder to adapter.
- Screw mains lamp to lampholder, connect it to power supply, and switch it on.

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


## 11.2 Low-voltage lamp 6 V/10 W

The 6 V/10 W low-voltage lamp is for the inclined incident illumination of spatial objects and is suitable for **observation and photomicrography**. The colour temperature of the halogen bulb is 2700 K.

For shadow-free illumination and for photomicrography, it is recommended that two low-voltage lamps be used; they can both be connected to the same 7 V regulating transformer.

The low-voltage lamp can usefully be connected to the **side-faced column**, so that during focusing the illumination is always directed on to the object (fig. 2).

The same advantage is gained if the low-voltage lamp and adapter (49) are connected to a **cast foot** (fig. 3) or to the transmitted-light base EB.

### Fitting to the side-faced column:

(fig. 2)

- Fit the adapter (53) either above or beneath the drive housing, depending on the working distance.
- Fit one or two clampable lampholders (46) to the adapter.
- Use the rotatable coupling (50) to attach the low-voltage lamp to the lampholder.
- Connect the low-voltage lamp to the 7 V regulating transformer (see section 14).
- Hold the lamp housing and slacken the clamping screw (47).
- Align the lamp to the object and retighten the clamping screw.
- By turning the lamp socket (51), focus the light spot as required.
- Two **filter holders** diam. 50 mm can be connected to each of the two pins on the lamp housing and swung in when needed. It is recommended that a **KGI heat-absorbing filter** be used to prevent overheating of the specimen.

To fit the low-voltage lamp to the objective mount, see section 11.1.

### Changing the bulb:

Note: Always use a cloth to hold the halogen bulb.

An arrow on the lamp mount and a line on the lamp housing indicate the position at which the two parts can be separated by pulling.

The lamp mount clicks automatically into position again when refitted.

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
### 11.3 Low-voltage lamp 6 V/20 W

The 6V/20W incident lamp is a powerful source of intense and very uniform illumination with a colour temperature of 3200 K. The **large lampholder** can be secured in any position. Hold the lamp housing when slackening the clamping knob (55).

#### Fitting:

- The **side-faced column** accepts two lamps by means of the adapter (53).
- The **cast foot** and the short posts on the **transmitted-light base EB** accept the lampholder either via the adapter (61) or directly with a thread connection.

#### Fitting the bulb:

- Remove the rear part of the lamp housing.
  - Hold the 6 V/20 W halogen bulb in a cloth and push it into position.
  - Close the housing with the wire clip.
- Connect the lamp with the 7 V regulating transformer (see section 14) and switch it on.

#### Centring the halogen bulb:

The centring must be repeated whenever the bulb is changed.




- Pull out the condenser lens (56).
- Cause the light spot to fall vertically on a white paper about 25 mm away.
- Using screw (57), focus the image of the filament and its mirror image (fig. 4a).
- Displace the filament images with the horizontal movement screw (58) so that they are opposite (fig. 4b).
- Using the vertical movement screw (59), bring the two images into contact (fig. 4c).

#### Focusing the light spot:

When the condenser lens is fully pulled out, a sharply-defined illuminated field of diameter 16 mm is formed about 75 mm in front of it. By retracting the lens, the diameter can be increased to about 40 mm without changing the position of the lamp housing. A frequently-used mean position is denoted by a click stop.

#### Filters:

A **heat-absorbing filter** is built into the lamp. A **diffusing filter in swing-out mount** (60) is attachable to one of the rods. A maximum of four 32 mm diameter filters can be used.

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
#### 11.4 Ring illuminator

The ring illuminator, which produces light of daylight quality, is used for the shadow-free illumination of high-contrast spatial objects. The intensity suffices for magnifications of up to 40 $\times$ .

The ring illuminator cannot be used on the drive housing "S".


##### Fitting

- Remove the objective.
- Slide the holder ring over the objective mount.
- Replace the main objective.
- Using the screw (42), secure the lamp housing at the desired height.
- Set the choke to the **available voltage** as follows:
  - Uncrew the fuse holder (41).
  - Pull out the cover plate of the voltage selector and replace it so that the appropriate voltage is visible in the cutout:
    - 110 V position for 100–120 V supply
    - 220 V position for 200–240 V supply
  - Tighten the fuse holder.
- Connect the power cable.
- Switch on the illuminator with the toggle switch (40).

Fuses:

for 110 V: 200 mA, 5 $\times$ 20 mm (370879)

for 220 V: 500 mA, 5 $\times$ 20 mm (167650)

#### 11.5 Near-vertical illuminator

The housing of the near-vertical illuminator includes a deviating prism, which directs light on to the object at an angle of about 5° from the vertical. Depressions are fully illuminated, and tools no longer throw disturbing shadows.

**Please note** that you can use only the special  $f=100$  mm, 150 mm, 175 mm and 200 mm objectives in conjunction with the near-vertical illuminator.

##### Fitting

- Remove the adapter ring from the housing of the near-vertical illuminator and screw it on to the optics carrier in place of the objective.
- Secure the incident-light housing to the adapter ring by means of the clamping screw.
- Screw the special objective on to the housing.
- Attach the fibre-optic light guide to the left- or right side of the housing:
  - Slacken the clamping screw, remove the plug, and clamp the fibre-optic light guide.

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


## 11.6 Coaxial incident illuminator

Coaxial illumination is available for flat specular subjects such as integrated circuits or metal sections. The light is directed along the two beam paths.

Built-in polarising filters eliminate stray reflections.

If the **drive housing "S"** is used in the **vertical-beam position**, a **rotatable filter holder** with quarter-wave plate must be fitted in front of the objective (fig. 2). When it is turned, the degree of extinction is changed. The filterholder is only usable with achromatic objectives and does not fit over the plano objective.

The coaxial incident lamp housing is available either with an incorporated 6 V/20 W halogen lamp or with a socket for a Volpi fibre-optic light guide. It forms a compact unit with the stereomicroscope, and fits between the optics carrier and the phototube (see section 21).

The use of the low trinocular tube (fig. 3) is recommended.

- Remove the binocular tube (see section 8).
- Fit the coaxial illuminator to the optics carrier so that the lamp housing is on the right and pointing backwards at 45°.

Manoeuvre it into position, ensuring that the locating screw engages the guide, and tighten it.

- Position the binocular tube on the coaxial illuminator and secure it.

For photography with the coaxial illuminator, a **quartz plate** must be added inside the camera.

**Detailed instructions** for the use of the coaxial illuminator are to be found in instruction booklet M2 232.

### Note:

The coaxial illuminator introduces a magnification factor of 1.5x.

Please note that large fields of view are not fully illuminated.

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


## 12. Transmitted light

The side-faced column is attached to the transmitted-light stands by means of an adapter (see section 7).

### 12.1 Transmitted-light stand, bright field, circular

The circular transmitted-light stand is suitable for the examination of transparent objects in bright field. For photography, it is recommended that the transmitted-light stand EB or the transmitted-light stand for bright and dark field be chosen.

In the circular transmitted-light stand for bright field, the light from the 6 V/10 W halogen bulb is directed through the specimen and into the objective by means of a deviating mirror, a collector lens, and a condenser lens. The colour temperature is 2700 K.

- Assemble the outfit as described in sections 7-10.
- The cable is packed underneath the stand; connect it to the 7 V regulating transformer (see section 14).
- Move the lever (30) to its limit and extract the guide (34) completely.

Hold the 6 V/10 W halogen bulb with a cloth, insert it, and replace the guide.

If even more light is required for special applications, a 6 V/20 W halogen bulb can be inserted.

The light spot is focused by displacing the guide with the lever (30) so that each field of view can be optimally illuminated. Since the illumination is very bright, it is however recommended that the entire glass stage plate (diam. 80 mm) be illuminated.

A diffusing filter is present in front of the collector lens. If required, a different 32 mm diam. filter can be fitted instead:

- Pull out the lamp cable.
- Loosen the three screws (33) on the underside of the base and remove the cover plate.
- Remove the spring clip from the filter-holder (32).
- Substitute the filter and replace the spring clip.
- Replace the cover plate.

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
### 12.2 Transmitted-light stand EB

A complete illumination system including a 6 V/20 W halogen bulb, a regulating transformer and a thermal fuse, is built into the transmitted-light stand EB (fig. 1).

The specially-designed collector system permits optimal use of the light. This stand is particularly useful for **photomicrography in transmitted-light bright field**. The colour temperature is 3200 K.

**Further details** about the operation of this stand and in particular about the use of the interchangeable collector lenses and the swing-in frosted filter are given in the instruction booklet M2 295.

### 12.3 Transmitted-light stand, bright/dark field

The transmitted-light stand for bright and dark field (fig. 2) enables **observations and photography to be carried out in either bright field or dark field**.


In the dark field technique, no light passes directly through the specimen into the objective. Light reflected at the interfaces of components with differing refractive indices travels into the objective, however, so that a bright image on a dark background is seen and the contrast of weakly-coloured or transparent specimens can be increased.

The built-in 12 V/100 W halogen bulb produces high light intensity at a colour temperature ideal for photomicrography (up to 3400 K).

It is connected to the 12 V regulating transformer, which simultaneously controls the fan.

The changeover from bright to dark field and back is carried out with the lever (35) at the rear of the base.

**Further details** on the use of this stand are to be found in the instruction booklet M2 280.

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

### 13. Combined illuminators

Numerous objects simultaneously display many of their structures only when they are illuminated by a combination of lighting techniques.

**Combined inclined and coaxial incident light** (fig. 1) enables objects possessing both flat specular areas and uneven parts, such as hybrid circuits and mechanical components, to be inspected.

- A 6 V/10 W or 6 V/20 W low-voltage lamp and the 6 V/20 W coaxial illuminator can be connected simultaneously to a 7 V regulating transformer (see section 14) which accepts a maximum loading of 50 W.

The **combination of incident and transmitted light** (fig. 2) is recommended for partly-transparent objects such as fibres, insects, and textiles. The inclined incident light emphasises details and structures of the surface.

- The 7 V regulating transformer for the circular transmitted-light stand will accept a low-voltage lamp in addition.
- For photomicrography, the low-voltage lamps should be connected to a separate transformer, in order that the colour temperature of the halogen lamp remains constant. The same is valid for the combination of the low-voltage lamp with the transmitted-light stand for bright and dark field.

The incorporated regulating transformer of the transmitted-light base EB will power the transmitted-light source and an incident lamp simultaneously.

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




## 14. Regulating transformer

The 0–7 V/50 VA **regulating transformer** (fig. 1) will accept the 6 V/10 W and 6 V/20 W low-voltage lamps, the circular stand for transmitted light with 6 V/10 W halogen bulb and the coaxial illuminator with 6 V/20 W halogen bulb. A maximum loading of 50 W is to be connected to the two sockets (66) (fig. 2).

The 0–12 V/100 VA **regulating transformer** (fig. 3) is used with the transmitted-light stand for bright and dark field with 12 V/100 W halogen bulb and with the coaxial illuminator if a 12 V/20 W halogen bulb is fitted. It has three sockets which accept a maximum combined loading of 100 W (fig. 4).

- Set the **voltage selector** (67) to the voltage available:
  - Position 115 for 100–120 V
  - Position 220 for 200–240 V
- Connect the power cable to socket (65) and to the power supply.

### Setting the limiting voltage:

- The click stop 6 V or 12 V prevents overloading of the bulb.
- Turn the regulating knob (63) clockwise until it engages.
- Slacken the screw on the metal ring (64).
- Turn the regulating knob back to zero.
- Connect the lamp to the socket (66). Turn the regulating knob until the voltmeter shows 6 V or 12 V, as appropriate.
- Tighten the screw on the metal ring.

### Changing the fuse:

- Pull out the power cable.
- Remove the screw (68) from the underside.
- Lift off the housing.
- Replace the defective fuse:
  - 7 V: 8 A, 5×20 (214477)
  - 12 V: 2×800 mA, 5×20 mm (343950)

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

## IV. Accessories

### 15. Filterholder

If the coaxial incident light housing is used in conjunction with the drive housing "S", a **filterholder with quarter-wave filter** must always be added (see section 11.5).

A **filter slide housing** is also available. This will fit between optics carrier and photo- or binocular tube (fig. 1).

The filter slide accepts two filters which can be readily changed. Ordinary gelatine filters can be cut into strips to match the filter hole size of 23 mm×45 mm.

### 16. Polarisation

Double-refracting substances, such as many crystals, petrological thin sections and fibres, can be more easily seen and identified in polarised light.

The transmitted-light stands accept either the **polariser on glass stage plate** or the **rotatable Pol. stage**. All polarising sets require an **analyser** in rotatable mount, which is fitted to the objective mount.

If photographs are taken in polarised light, a **quartz plate** must be added inside the shutter-piece.

#### 16.1 Polariser on glass stage plate

- Place the polariser on glass stage plate, polariser downwards, in the cutout of the transmitted-light stand (fig. 3).
- Rotate the analyser using the lever (83) so that the field of view without specimen is as dark as possible.
- Place the specimen in position and focus on it.

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


## 16.2 Rotatable Pol. stage

- 83 lever for rotatable analyser
- 84 centring screws
- 85 lever to swing in polariser
- 86 screw to secure rotatable stage
- 87 milled ring to rotate compensator

The rotatable Pol. stage is equipped with a swing-out **polariser**.

If required, a **first-order red compensator** (sensitive-tint plate) and a **mechanical stage** can be added to it. The compensator facilitates the identification of weakly-birefringent specimens.

- To fit the compensator, remove the screw (88) from the underside of the stage.
- Push the compensator in its metal mount fully into the stage.
- Retighten the screw.
- Fit the mechanical stage using the screw (82).
- Remove the glass stage plate and the stage clips from the transmitted-light stand.
- Fit the rotatable Pol. stage so that the two centring screws (84) are symmetrically positioned in relation to the column.
- Secure the stage with the lever (81).
- Switch on the light.
- Swing in the polariser with the lever (85) and pull out the compensator slide.
- Look into the binocular tube and move the lever (83) until the field of view without specimen is dark as possible.
- Secure the specimen with the arm of the mechanical stage and focus on it.

With the help of a crosshair in the eyepiece (see section 10) and of the two centring screws, **centre** the stage as follows:

- The point about which the stage rotates is adjustable and must be brought into coincidence with the axis of the objective, i.e. with the centre of the crosshair.
- Spin the stage in order to establish the point of rotation in the field of view.
- Using the screws (84), displace the point of rotation in the directions ab and cd (fig. 3). Successively repeat this procedure, simultaneously rotating the stage, until the point of rotation coincides with the centre of the crosshair.
- Push in the compensator mount and rotate it until the field of view appears red. The rotatable Pol. stage is calibrated and can be used for **angular measurements** (see instruction booklet M2275).

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




## 17. Stages

The use of special stages simplifies work with Wild stereomicroscopes.

They are attachable to all incident- and transmitted-light stages with 80 mm diameter port, to the transmitted-light stage, and to stage carriers with magnetic linkage (fig. 2). Insert the spacing ring (89) if the trapezoidal incident-light stand is used (fig. 1).

- Remove the stage plate and stage clips supplied as standard (fig. 1).
- Reposition the stereomicroscope correspondingly higher on the side-faced column.

### 17.1 Gliding stage

The gliding stage (fig. 3) can be quickly moved in any direction and rotated.

- Fit the stage and secure it with the lever (81).
- Fit a glass or black/white stage plate and the two stage clips.

### 17.2 Cup stage

The cup stage (fig. 4) is essential for examining cultures in petri dishes which are held in position by a movable holder.

Insects and plants are secured to the rubber surface with needles.

- Fit the supporting ring.
- Add the cup stage, which can be tilted to any position.

The freedom of movement of the cup stage can be further increased by mounting it on the gliding stage.

### 17.3 Stage carrier with mechanical stage

The stage carrier with mechanical stage (fig. 5) facilitates the systematic scanning of specimens.

- Fit the stage carrier and secure it with the lever (81).
- Add the appropriate glass stage plate.
- Secure the mechanical stage to the stage carrier with the Allen screw (82). The mechanical stage can also be used on the rotatable Pol. stage.

The object is secured with the arm and is moved in an x- or a y-direction with the knobs of the mechanical stage (76 mm x 28 mm).

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

#### 17.4 Transmitted-light stage

Place the spacing ring and the transmitted-light stage in the stage cutout (fig. 1). Fit a glass (or acrylic) stage plate and the two stage clips. The light source can be daylight or the 6 V/10 W low-voltage lamp. The lamp-holder can be attached to the post. Cause the light spot to fall on the tiltable mirror. Slight deliberate disadjustment of the mirror will produce inclined transmitted illumination, which is very advantageous for observing small translucent objects.


#### 18. Handrests

The individually-adjustable handrests ensure precise and fatigue-free work on the object (fig. 2). A rubber ring on the base prevents sliding on the table surface. The upper surface of the handrest is insulated with a layer of paint.

#### 19. Attachment for vertical and oblique observation

This attachment enables three-dimensional objects to be viewed and photographed either vertically or from all around at an angle of 45° without it being necessary to turn or tilt them (fig. 3). The attachment can be attached to the mount of the 1.0× achromatic objective by means of the securing ring (13).

Further details are given in the instruction booklet M2 265.

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
## 20. Measuring

Because of its high optical quality, the Stereomicroscopes of the WILD M3 Series are particularly suitable for precise measurements.

For measurement work, adjustable eyepieces can be furnished with a **graticule (reticle) in mount** (see section 10). Both eyepieces must then be of the adjustable type (fig. 1).

When using the drive housing "S", position the measuring eyepiece in the beam path which is to be used for measurement.

The graticules available have various scales or grids, or a crosshair.

Angles can be measured with the **10× wide-field goniometer eyepiece**, or with the **rotatable stage** used for polarisation work.

**Further details** of procedures used for measuring are given in the instruction booklet M2 275.

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
## 21. Tubes

The stereomicroscopes of the WILD M3 Series have a removable binocular tube, so that various accessories such as phototubes, the beam splitter, the drawing tube or the double-iris diaphragm can be fitted.

- Remove the binocular tube (see section 8).
- Mount the accessory tube.
- Replace the binocular tube on the accessory tube.

It is advantageous to use the **low binocular tube** with tall combinations.

### 21.1 Double-iris diaphragm

With the double-iris diaphragm, which fits between the optics carrier and the binocular tube (fig. 1) or the phototube A, the depth of field can be matched to the characteristics of the object (see section 22.5). This is particularly advantageous in photomicrography in monocular assembly, because no double-iris diaphragm is incorporated in the phototube A.

Double-iris diaphragms are however an integral part of the **binocular tubes T**, of the **trinocular phototubes** and of the **beam splitters**.

When the diaphragms are partly closed, the light cone entering the objective is reduced after passing through the magnification changer. The depth of field increases and the image brightness is reduced.

The iris diaphragms are operated with the lever (71).


Position 1 is the smallest aperture/position; 'OPEN' is the largest.

### 21.2 Discussion tube

The stereomicroscopes of the WILD M3 Series can be expanded into excellent training instruments by combining them with a discussion tube (fig. 2). Two observers see the same image simultaneously, so that both instruction and discussion are possible.

- The discussion tube is connected to the universal stand by means of a drive housing.
- The optics carrier is fitted beneath the discussion tube, and the two binocular tubes are fitted above it.
- The objective adapter (75), which also serves as a lampholder, is connected to the optics carrier by means of two hollow screws. The low-voltage lamps are secured using an adapter diam. 25 mm.

**Further details** are given in the instruction booklet M2 261.

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
### 21.3 Drawing tube

In a drawing, characteristic features can be emphasised, less important detail can be omitted, and features in various planes of the object can be depicted. The drawing tube can be used in either incident or transmitted-light specimens in daylight and in a comfortable sitting position. The user sees the object and the drawing simultaneously. The mirror is protected from damage by a quartz coating.

#### Fitting

- Fit the drawing tube between the optics carrier and the binocular tube so that the mirror points to the left or to the right.
- Move the lever (72) to the white spot so that the drawing surface is imaged in the binocular.
- Move the objective (74) of the drawing tube so that both the drawing surface and the pencil are seen in focus. The focusing range is from 100 mm to infinity. It is recommended that the drawing surface be brightened by means of a separate table lamp, e. g. with a Wild 25 W mains lamp on cast foot.
- Place the object in position, switch on the illumination, and bring the image into focus.
- The red screw (73) releases the mirror which can be removed and turned in order to mirror in the images of documents or objects in various positions. Comparisons can thus easily be made.

If individually-made scales or grids are mirrored in and photographed with the object, measurements of number, length and area can be carried out.

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
## 21.4 Phototubes

**Phototube A** is for a monocular camera assembly; **Phototubes HU and HV**, and the **inclined trinocular tube**, are for a trinocular camera assembly. The **trinocular tube** is a combined tube for observation and photography (fig. 1). It permits any camera to be attached. The low viewing height reduces the overall height of the outfit.

The HV Phototube has a beam splitter system which can direct all of the light to the binocular or can deflect half of it to the camera. The HU Phototube can direct all of the light to the binocular or can direct to the camera all of the light reaching one of the eyepiece tubes of the binocular.

- Position the phototube over the left or the right beam path of the optics carrier and secure it.

When the drive housing "S" is used, the phototube is to be positioned over the beam path which will be used for photography.

- Accessories such as the coaxial incident illuminator, the double-iris diaphragm and the filter slide housing are to be fitted between the optics carrier and the phototube.

- Equip the HU and HV phototubes with a binocular tube and with two 10 $\times$ , 15 $\times$  or 20 $\times$  adjustable eyepieces.

One eyepiece requires a format-indicating graticule.

### Camera for Wild Modular Photomicrographic System

To fit a Wild shutterpiece, first place an eyepiece tube with a fixed eyepiece in the phototube.

**Detailed information** is available in the instructions for your Wild Modular Photomicrographic System.

### Single-lens reflex camera

Single-lens reflex cameras can be connected to Wild stereomicroscopes via a phototube with eyepiece tube, fixed eyepiece, connecting tube, 0.32 $\times$  camera objective, adapter and T-mount.

### TV and cine cameras

For **single-stage imaging**, the connecting member is an adapter with C-mount.

For **two-stage imaging**, the connecting member consist of an eyepiece tube with fixed eyepiece and an objective  $f = 50$  mm with C-mount.

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
### 21.5 Beam splitter

Two versions of the beam splitter, with different beam partition characteristics, are available:

- for dual observation and television (**50% to the binocular** and 50% laterally).
- for photo and cine (**20% to the binocular** and 80% laterally).

The double-iris diaphragm, which is built in, influences the beams in the binocular and in both of the lateral outlets. The latter can be used simultaneously for different purposes (fig. 1).

- Fit the beam splitter between optics carrier and binocular tube so that the diaphragm lever (71) points towards the observer.
- The phototube, the cine/TV tube and the tubes for a second observer have a self-centring connection to the dovetail mount of the beam splitter.

### Phototube

The phototube (fig. 2) can be used with commercially-available single-lens reflex cameras. A T-mount (100) attaches these to the  $f = 250$  mm or  $f = 350$  mm adapters.

### Cine/TV tube

The cine/TV tube (fig. 3) can be used for TV and cine cameras which accept a C-mount (101). The four possible engagement positions of the cine tube enable the camera to be positioned where it is least in the way.

A **format-indicating graticule** (reticle) for the phototube and for the cine/TV tube indicates the format limits for 35 mm film for adapter tube lengths of 250 mm and 350 mm respectively. The graticule fits into the adjustable wide-field eyepieces and into the 10 $\times$  eyepieces for spectacle wearers (see section 10). If a graticule is used, **two adjustable eyepieces** are required.



### Tube for second observer

This tube is available in a **short** (fig. 1) and in a **long version**.

Because the tube can be rotated about two axes, it can be positioned at any angle. The knurled sleeve (104) enables the position of the image to be corrected.

### Stereo attachment for second observer

The advantage of this accessory is that the tube can be moved to any position and that the second observer enjoys comfortable binocular observation (see page 72).

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## V. Appendix

### 22. Choosing the optics combination

The choice of optics combination is influenced by the **working distance** required, the **total magnification** obtainable, and the **appropriate field of view**.

**Resolution and depth of field** are also dependent on the optical components chosen.

The table (page 79) gives field diameters and total magnifications for the three models of the M3 Series when used in conjunction with the wide-field eyepieces 10×, 15×, 20× and 32× and the objectives 1.0×, 0.32×, 0.5×, 0.63×, 0.8×, 1.5× and 2.0×.

#### 22.1 Working distance

The working distance can only be changed by substituting or adding an **objective** of different focal length.

The higher the objective magnification, the shorter the focal length and the working distance.

<u>Achromatic objectives</u>	<u>Working distance</u>
0.32×	297mm
0.5×	187mm
0.63×	149mm
0.8×	112mm
1.0×	89mm
1.5×	49mm
2.0×	27mm

#### Planachromat objectives

1.0×	89mm/104mm*
1.6×	28mm

\* without adapter ring for accessories

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
## 22.2 Total magnification

The imaging within a stereomicroscope takes place in two stages.

In the **first stage**, the objective together with the magnification changer forms the intermediate image of the object.

In the **second stage**, the intermediate image is further magnified by the eyepiece.

The **total magnification is the product of the magnifications of eyepiece and objective**, and is calculated as follows:

$$\frac{M_{ob} \times M_{oc} \times p \times t}{10} = M_{total}$$

Example:  $\frac{1.5 \times 15 \times 16 \times 1.5}{10} = 54 \times$

- $M_{ob}$  magnification of objective
- $M_{oc}$  magnification of eyepiece
- $p$  magnification changer step
- $t$  tube factor (e.g. 1.5× for coaxial illuminator)
- $M_{total}$  total magnification

The magnification range can be extended by using various additional objectives. The 0.32× and the 2.0× objectives produce lower and higher ranges respectively.

The 10×, 15×, 20× and 32× eyepieces also influence the magnification, so that total magnifications from 1.9× to 256× are attainable (see table p. 79).

## 22.3 Field diameter

The field of view is the circular area of the specimen which can be seen down the microscope (fig. 1).

It is **influenced by objective and eyepiece characteristics** and decreases with increasing magnification, e.g. 111 mm at 1.9× and 1.5 mm at 160× (see table p. 79).

The field diameter is calculated as follows:

$$\varnothing = \frac{\text{field number} \times 10}{M_{ob} \times p \times t}$$

The field numbers are inscribed on the eyepieces (e.g. 10×/21; field number: 21).

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


## 22.4 Resolution

The resolving power of a stereomicroscope is a measure of its **ability to distinguish fine detail**, i.e. to separate two closely-spaced points.

Resolution and contrast depend on the **objective and on the magnification changer position**. Resolution is proportional to the angle of aperture (fig. 1). Powerful objectives of short working distance thus give high resolution, and this increases with the zoom position.

**Details not resolved by the objective cannot subsequently be rendered visible by using strong eyepieces.**

A given **total magnification** is attainable either by using a weak objective coupled to a strong eyepiece, or vice versa:

e.g. 80× with a 1.0× objective, 20× eyepieces and step position 40, or with a 2.0× objective, 10× eyepieces and step position 40.

The second combination is preferable if the working distance is unimportant.

An additional eyepiece magnification should be so chosen that all details resolved by the objective are recognisable. No advantage is gained by an excessive secondary magnification, but an inadequate secondary magnification means that the potential of the objective is not being fully exploited.

## 22.5 Depth of field

Depth of field depends on the objectives, the eyepieces and the tube factor.

**Depth of field decreases with increasing magnification.**

The user can use these laws to match the appropriate optics combination to the characteristics of the specimen.

For larger objects with coarse structures, the double-iris diaphragm can often be used to improve depth of field at the expense of resolution (see section 21.1).

Fig. 2: Double-iris diaphragm open

Fig. 3: Double-iris diaphragm partially closed

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
*WILD M3Z Stereomicroscope  
with stereo attachment for second observer*

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### 23. Care of instrument

Wild Stereomicroscopes are high-performance precision instruments of the highest optical and mechanical quality. The security offered is correspondingly high.

**We guarantee the quality of each of our instruments. Our guarantee covers all faults in materials and manufacture. It does not, however, cover damage resulting from careless or improper handling.**

It is evident that a valuable optical instrument merits considerate treatment. It will then fulfil its functions for decades with the same precision.

**Dust and dirt** have a very bad effect on the image quality.

- Use the dust cover therefore, during pauses in work.
- Eyepieces, tubes without eyepieces, and the optics carrier when without tube, should always be protected against dust by means of the covers provided.
- **Dust on optical components** is best removed with a blow-bulb or with a soft brush.
- Tenacious marks on objectives and eyepieces can be removed with a soft lint-free cloth dipped in pure alcohol. Special cloths for cleaning optics are also available.
- Never dismantle optical systems unless instructions for doing so are given in the instruction booklet.
- Protect optical instruments from **damp, fumes, corrosive substances, and acid or alkaline vapours**. Do not store chemicals nearby.
- The **mechanical parts** of the instrument can be cleaned with a soft brush or with a soft cloth.
- Never grease, oil or dismantle mechanical parts unless instructions for doing so are given in the instruction booklet.

If the instrument no longer functions satisfactorily, call in a specialist or your local Wild Leitz agent, or contact the factory.

Modifications resulting from technical developments may be made in the interest of our customers. Thus, illustrations, specifications and catalogue references are not binding and are subject to change without notice.

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


\*-89 mm with 1.0× achromatic objectives and 1.0× planachromatic objective with adapter ring for accessories.  
-104 mm with 1.0× planachromatic objective without adapter ring.

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Illustrations, descriptions, and technical data are  
not binding and may be changed without notice.

M2 131 III 90 - Printed in Switzerland  
© Wild Leitz Ltd, Heerbrugg

**WILD**  
LEITZ

Trademark of world-famous products  
from the Wild Leitz Group